

Expander

- The expander is an orthodontic appliance that is cemented to the teeth. It is not taken in and out by the patient.
- The expander spans across the roof the mouth. The middle portion of the appliance has a screw in it. When turned, the expander gets wider.
- When the expander is cemented, we will instruct you on how to turn the screw and thereby activate the expander.
- Steps in making an expander:

Appointment #1 – Place separators between the teeth (10 minutes)

Appointment #2 – One week after separators, fit stainless steel bands around the molars and take impression of the teeth (45 minutes) You will leave with new spacers.

Appointment #3 – One week after fitting the bands, cement the expander in place (45 minutes)

- The expander is checked every 2 weeks to determine how many turns will be needed.
- Once there have been enough turns, the appliance will remain in the mouth for four months to hold the expansion.



PALATAL EXPANDER

Your child has been fitted with a special appliance made to widen the upper arch. Please follow the instructions below. If you should have any questions, please feel free to contact our office.

- 1. Use the key to turn the appliance as prescribed. Your assistant will make sure to explain how many turns to complete between appointments.
- 2. Insert the key into the front hole and <u>gently</u> push steadily toward the back of the throat. This will bring the next hole into view. **DO NOT BACK UP**. Remove the key toward the throat.
- 3. The patient will feel pressure for a few hours each time the key is turned during the first few days. Take Tylenol, or other similar products, if needed.
- 4. A space will open between the two front teeth within a few days. This is normal and nothing to cause concern. It will close by itself over time.
- 5. If you forget a turn, **DO NOT** give extra turns to try to catch up. Just continue with the prescribed schedule.
- 6. DO NOT EAT HARD, CRUNCHY, OR STICKY FOODS.
- 7. Keep the appliance clean. The patient should brush around all parts of the appliance, especially next to the gums.
- 8. If one of the bands comes loose, or you have any questions, please call our office.

